

# SaeboMAS

The **SaeboMAS** dynamic mobile arm support system is a zero gravity upper extremity device specifically designed to facilitate and challenge, the weakened shoulder and elbow during functional tasks and exercise drills.

Patients suffering from neurological and orthopedic conditions and who exhibit proximal weakness can benefit from the much needed support and facilitation. Performing exercise drills with minimal compensation, as well as self-care tasks like eating or drinking, are just a few of the activities patients will be able to achieve with the **SaeboMAS**.

Additionally, patients with proximal weakness can now participate in proven treatment techniques of highly repetitive, task oriented activities which would have otherwise been impossible.



## Features:

- ▲ Adjustable spring based parallelogram offering various levels of assistance.
- ▲ Incorporate all planes of movement so multi-directional activities can be achieved.
- ▲ Measurable graded Tension Scale for tracking and documenting progress.
- ▲ Height adjustable.
- ▲ Table Mount.
- ▲ Comfortable malleable forearm support with removable liners for infection control.
- ▲ Elbow Support attachment.
- ▲ Lightweight and portable.



## Benefits:

- ▲ Increase motor control, strength and ROM.
- ▲ Improve ADL performance.
- ▲ Reduce hand-over-hand assistance, allowing the patient to take on a more active role.
- ▲ Minimize/prevent over-use injuries and unwanted movement.
- ▲ Provide the clinician with an “extra pair of hands” to offer more effective facilitation when needed.
- ▲ Provide an opportunity for the patient to perform highly repetitive task-oriented drills.
- ▲ Safe and effective way to treat shoulder subluxation.

## Combining SaeboMAS with SaeboFLEX<sup>®</sup>

