

# SaeboStretch<sup>®</sup>



Finally, a **Dynamic Solution**  
for a **Dynamic Problem**<sup>™</sup>

## What is the SaeboStretch?

## How is the SaeboStretch different from current splints?



The SaeboStretch® is a revolutionary new award winning dynamic hand splint for individuals suffering from neurological injuries such as stroke. The splint helps to prevent joint damage while improving and maintaining range of motion.

### Positions the Hand:

The SaeboStretch provides maximum positioning utilizing the following:

- Proprietary new dynamic hand piece that comes in three varying grades of resistance to better match your patient's specific needs.
- Proprietary new strapping system that utilizes key anatomical points of control to maximize control and insure an intimate fit.
- Malleable wrist and thumb mount allows the therapist to further customize the fit.
- New non-slip strap material and design helps facilitate keeping the fingers in the desired position.
- Proprietary Breeth-O-Preen covering that allows moisture to wick away enhancing skin integrity and hygiene.
- Special padding to help maintain the integrity of the palmer arch.




other splint



SaeboStretch

### Protects the Joints:

Traditional splints are made from hard/rigid material that do not allow for any movement of the fingers. The SaeboStretch has a unique dynamic hand piece that allows the fingers to move through flexion caused by postural changes, associated reactions and increased tone. The same stretch technology then gradually repositions the fingers into extension. The SaeboStretch includes three different hand pieces that offer various grades of resistance so you can further match your patient's specific needs.



*"The SaeboStretch finally allows the control of my patient's hand that I have always wanted. He stood, he walked, and even yawned without losing position of his fingers. Thanks for making the product all therapists have desired."*

**Liz Miner**  
Occupational Therapist

*"It really is the most comfortable hand splint I've ever worn. Had the SaeboStretch been available when I had my stroke, I think I might have more function and use of my hand today."*

**John Whitfield**  
Stroke Survivor

## What happens when you put a dynamic hand in a static splint?

- Joint Deformities
- Contractures
- Hypermobility



Give your patient's hand the best possible chance for recovery

When choosing a splint, request the best. Request the **SaeboStretch**.

To learn more about this revolutionary new splint, visit our Website at [www.saebo.com](http://www.saebo.com) or call **888-284-5433**.

  
**Saebo**  
Leaders in Stroke Rehabilitation

With **Saebo**  
There's No Plateau In Sight.®



2725 Water Ridge Parkway  
Suite 320  
Charlotte, NC 28217  
888.284.5433  
[www.saebo.com](http://www.saebo.com)