1. Find Your "Splint" Size

(A) Wrist Circumference: measure around wrist just below wrist bone.





(B) DORSAL HAND WIDTH

(B) Hand Width: measure dorsal hand width just below knuckles of your pointer and pinky fingers (do not wrap tape down side of hand); use a straight ruler if possible.

			Wrist Splir	nt Sizing Ch	art							
	Wrist Circumference (inches)											
Dorsal Hand Width (inches)		4 - 4.5	4.5 - 5.5	5.5 - 6.5	6.5 - 7.5	7.8 - 8.5	> 8.5					
	< 2.5	S	S	S								
	2.5 - 3.0	s	S	М	М							
	3.0 - 3.5	S	М	М	М	М	L					
	3.5 - 4.0			М	М	L	L					
	> 4.0				L	L	L					

2. Find Your "Glove" Size

(A) Hand Circumference: measure around the hand, just below the knuckles.



(B) 3RD FINGER LENGTH

(B) 3rd Finger Length: measure from the tip of the middle finger down to the edge of the palm.



(A) HAND CIRCUMFERENCE



			Glove Lin	er Sizing C	hart							
	Hand Circumference (inches)											
hes)		4.5 - 5.5	5.5 - 6.5	6.5 - 7.5	7.5 - 8.5	8.5 - 9.5	9.5 -10.5					
3rd Finger Length (inches)	6.0 - 6.5	XS	XS	XS								
	6.5 - 7.0	XS	XS	XS	s							
	7.0 - 7.5		s	S	S	М						
3rd F	7.5 - 8.0			М	М	L	XL					
	> 8.0				L	XL	XL					