

Meet the New & Improved SaeboStretch



Dear Valued Customer,

Saebo is pleased to announce that we have made some exciting improvements to one of our most popular products. Thank you for your feedback. Your input has been essential in making this revolutionary splint even better.

What are the changes?

Straps:

- Upgraded straps now consist of comfortable elastic woven non-slip material.
- Straps are sewn to the cover to maximize proper positioning and ease of donning.

Cover:

- New cover features a zipper closure that can easily be removed for cleaning and changing hand plates.
- Cover also includes non-slip gel strategically located on the hand and thumb to minimize migration.

Thumb System:

- The new thumb system now offers 2 adjustments to further customize the thumb position. Position the thumb in palmar abduction/adduction and/or radial abduction/adduction.

Hand Plate Screws:

- The small hex screws have been replaced with standard size head cross screws (i.e., Phillips head) for improved performance.



*Dynamic Solution for
a Dynamic Problem™*