

Now Available in 5 Sizes

SM/MD/LG Sizes



- Malleable wrist and thumb supports
- Removable non-slip cover
- 3 energy-storing hand plates to adjust resistance
- Adjustable thumb system
- 2 forearm straps
- 1 hand strap
- 3 finger straps
- 1 thumb strap

Colors: black with blue accents

XS Size



- Malleable wrist and thumb supports
- Removable non-slip cover
- 3 energy-storing hand plates to adjust resistance
- Adjustable thumb system
- 2 forearm straps
- 1 hand strap
- 1 finger strap
- 1 thumb strap

*Colors: black with blue accents and blue trim
black with pink accents and pink trim*

XXS Size



- Malleable wrist and thumb supports
- Removable non-slip cover
- 1 energy-storing hand plate
- 2 forearm straps
- 1 hand strap
- 1 finger strap
- 1 thumb strap

*Colors: black with blue accents and blue trim
black with pink accents and pink trim*

**Adjustable thumb system not included*