

SaeboStim Pro



Clinically effective solution for muscle rehabilitation.

The **SaeboStim Pro** is an advanced Neuromuscular Electrical Stimulation System. The **SaeboStim Pro** was developed to prevent or slow down muscle weakness, improve motor impairment, and re-educate and strengthen weakened muscles. The device works by triggering muscles to contract and relax using electrical impulses. This portable, digital, dual channel stimulation device incorporates a comprehensive range of programs to aid in the rehabilitation process at home.

Features:

- Unique Trigger Mode for functional training. Stimulation triggered at any point during functional task by pressing Trigger Key.
- When Key is released, stimulation turns off.
- Large clear digital display with easy to understand symbols.
- One touch operation.
- Display screen shows length of time left/elapsed in the current session.
- 14 easy to use pre-set programs created to treat all levels of muscle weakness and many specific conditions.
- Sensory Electrical Stimulation (SES) program included.
- Full customization mode allowing for adjustable pulse rates, durations, ramp up/ramp down, on/off time controls.
- Consider use of **SaeboStim Pro** for conditions such as:
 - Disuse atrophy
 - Muscle spasms
 - Range of motion
 - Muscle re-education
 - Lack of blood circulation
 - Chronic or post-surgical pain



Simple, powerful, muscle stimulation made easy

1. On/Off Button
2. Program Selection
3. Lock Button
4. Trigger Button
5. Intensity Controls
6. Belt Clip (on back)
7. 2 Channel Output
8. LCD Screen