

SaeboStim Spa

Aquatic Electrical Stimulation Therapy



What is the SaeboStim Spa?

The **SaeboStim Spa** provides extremely safe, low-level electrical stimulation designed to reduce spasticity, improve sensation, minimize muscle weakness, improve function, re-educate and strengthen weakened muscles, and decrease pain. The device works by submerging the involved foot or hand, along with the associated rubber electrodes, into the water for a 30-minute stimulation session.

The **SaeboStim Spa** provides circumferential “under water” nerve stimulation that is consistent, comfortable and evenly distributed throughout the affected area. This form of stimulation surrounds the targeted tissue with a soothing and uniform level of current that treats the underlying conditions.

One Device. Multiple Conditions.

In addition to providing low-level sensory electrical stimulation for clients suffering from stroke and other neurological injuries, the **SaeboStim Spa's** specialized waveforms provides treatment for additional conditions including peripheral neuropathy and nerve pain disorders. This breakthrough device can calm the nerves and increase blood flow which allows the nerve cells the potential to recover.

